Instructions for Use

Instructions for Use - GB
Gebrauchsanweisung - DE
Instructions d'utilisation - FR
Istruzioni per l'uso - IT
Manual de instrucciones - ES
Gebruiksaanwijzing - NL
Instrukcja obsługi - PL

Impedance Audiometer Titan

80655311 ver. 06/2011

Instructions for Use - GB

Impedance Audiometer Titan


1.1 Intended Use

The Titan with IMP440 Impedance System is an electroacoustic test instrument that produces controlled levels of test tones and signals intended for use in conduction diagnostic hearing evaluations and assisting in the diagnosis of possible otologic disorders. It features tympanometry and acoustics reflex. The target population for Titan with IMP440 includes all ages.

The Titan with DPOAE440 is intended for use in the audiologic evaluation and documentation of ear disorders using Distortion Product Otoacoustic Emissions. The target population for Titan with DPOAE440 includes all ages.

The Titan with ABRIS440 is intended for use in the audiologic evaluation and documentation of ear and nerve disorders using auditory evoked potentials from the inner ear, the auditory nerve and the brainstem. The target population for Titan with ABRIS440 is newborns.

Typical Applications

Titan is well suited for screening, as well as diagnostic work. Newborn screening programs will particularly appreciate the Automated Auditory Brainstem Response (AABR) test implemented under the ABRIS module and the option to create user defined PASS/REFER protocols for DPOAE screening. The presence of the 1000Hz probe tone for tympanometry also allows more reliable tympanometry results in infants.

Titan, as a screener, will typically be operated as a hand held device. Results can be printed using the (Bluetooth) thermal printer or a Bluetooth PC printer. Optionally, measurements can be transferred later to a database on your PC for printing.

In combination with the cradle and shoulder box for IMP440 and DPOAE440, the Titan gives you more flexibility in placing the probe. Choice of ear, start and stop of the measurement can be controlled from the shoulder box.

In combination with the Preamplifier (PreAmp), the Titan is capable of performing the full range of tests from the ABRIS440, DPOAE440 and IMP440 modules.

The Titan, when licensed with a diagnostic or clinical version of the IMP440 or DPOAE440 software, is generally operated via the PC. PC operation allows access to advanced testing parameters, for example, it permits manual pressure control within the impedance module and the ability to add extra test points in DPOAE after testing. Testing is as quick in PC controlled mode as when using the handheld device. Results are saved in your database and can be printed on an external printer. Protocols can be modified with all the flexibility you need.

1.2 Safety Notes

 WARNING	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
 CAUTION	CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
 NOTICE	NOTICE is used to address practices not related to personal injury.

It is recommended that all staff working with the Titan read this manual and thereby learn the precautionary measures.

WARNING

1. External equipment intended for connection to signal input, signal output or other connector, shall comply with the relevant IEC standard (e.g. IEC 60950 for IT equipment and the IEC 60601 series for medical electrical equipment). In addition, all such combinations - *systems* - shall comply with the standard 60601-1-1, *Safety requirements for medical electrical systems*. Equipment not complying with IEC 60601 shall be kept outside the patient environment, as defined in the standard (at least 1.5 m from the patient).
2. Any person who connects external equipment to signal input, signal output or other connectors has created a system and is therefore responsible for the system complying with the requirements of IEC 60601-1-1. If in doubt, contact your service technician or local representative for help.
3. In order to comply with IEC 60601, use of the direct power adapter to mains for the thermal printer when it is connected to the cradle is not permitted. The cradle connection provides power for the thermal printer.
4. In order to maintain a high level of safety it is necessary to have the instrument and its power supply checked according to the medical electrical safety standard IEC 60601-1 on a yearly basis by a qualified service technician.

CAUTION

1. Use this device only as described in this manual.
2. Use only the disposable Interacoustics ear tips designed for use with this instrument.
3. Never insert the probe tip into the ear canal without affixing an ear tip as omission may damage the patient's ear canal.
4. Keep the box of ear tips outside the reach of the patient.
5. Be sure to insert the probe tip in a way which will assure an airtight fit without causing any harm to the patient. Use of a correct and clean ear tip is mandatory.

6. Be sure to use only stimulation intensities acceptable to the patient.
7. When presenting contralateral stimuli using the insert phones – do not insert the phones or in any way try to conduct measurements, without a correct insert ear tip in place.
8. When presenting contralateral stimuli using headphones – do not try to conduct measurements without the MX41 cushion in place.
9. Clean the headphone cushion regularly using a recognized disinfectant.
10. We recommend using a new ear tip for each patient.

NOTICE

1. Careful handling of the instrument whenever in contact with a patient should be given high priority. Calm and stable positioning while testing is preferred for optimal accuracy.
2. The Titan should be operated in a quiet environment, so that measurements are not influenced by outside acoustic noises. This may be determined by an appropriately skilled person trained in acoustics. OSHA defines a quiet room for audiometric hearing testing in its guideline OSHA 29 CFR 1910.95.
3. It is recommended that the instrument be operated within an ambient temperature range of 15°C / 59°F – 35°C / 95°F.
4. Storage in temperatures below 0°C / 32°F and above 50°C / 122°F may cause permanent damage to the instrument and its accessories.
5. The headphone and insert phone are calibrated to the shoulder box of this instrument – introducing transducers from other equipment requires a re-calibration.
6. Never clean the transducer housing with water or insert non-specified instruments into the transducer.
7. Do not drop and avoid other undue impact to this device. If the instrument is dropped or otherwise damaged, return it to the manufacturer for repair and/or calibration. Do not use the instrument if any damage is suspected.
8. Although the instrument fulfils the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment, caution must be taken to observe that no mutual disturbance appears.
9. The temperature of the battery pack in the thermal printer will rise if it is powered through both the printer's power supply and through the RJ12 cable to the cradle. Use either one or the other power connections, but never both, as this may severely reduce the lifetime of the battery.
10. The connectors marked with the symbol shown below should not be touched without taking necessary ESD (electrostatic discharge) precautions.
11. The best way to protect the Titan from ESD is to avoid the build-up of electrostatic charges through the use of conductive footwear and floor coverings. Humidification of the air is also an effective preventive means. Failing to comply with these measures can cause the build-up of voltage levels of several thousand volts. If the Titan is subjected to such an ESD pulse, it may be necessary to remove the battery for a few seconds and then restart the instrument.

The precautionary symbol:

Within the European Union it is illegal to dispose of electrical and electronic waste as unsorted municipal waste. Electrical and electronic waste may contain hazardous substances and therefore have to be disposed of separately. Such products will be marked with the crossed-out wheelie-bin image shown to the left. User cooperation is important in order to ensure a high level of reuse and recycling of electrical and electronic waste. Failure to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

1.3 General Maintenance Procedures

The performance and safety of the instrument will be maintained if the following recommendations for care and maintenance are observed:

1. It is recommended that the instrument go through at least one annual service, to ensure that the acoustical, electrical and mechanical properties are correct. This should be carried out by an authorized repairer in order to guarantee proper service and repair.
2. Before connection to the mains network, be sure that the local mains voltage corresponds to the voltage labeled on the instrument. Always disconnect the power cord if the instrument is to be opened for inspection.
3. Observe that no damage is present to the insulation of the mains cable or the connectors and that it is not exposed to any kind of mechanical load that could involve damage.
4. For maximum electrical safety, turn off the power from a mains powered instrument when it is left unused.
5. Do not place the instrument next to a heat source of any kind, and allow sufficient space around the instrument to ensure proper ventilation.
6. To ensure that the reliability of the instrument is maintained, we recommend that the operator at short intervals, for instance once a day, performs a test on a person with known data. This person could be the operator.
7. If the surface of the instrument or parts of it are contaminated, it can be cleaned using a soft cloth moistened with a mild solution of water and detergent or similar. The use of organic solvents and aromatic oils must be avoided. Always disconnect the mains power adaptor and battery during the cleaning process, and be careful that no fluid enters the inside of the instrument or accessories.
8. After each patient examination, ensure that there has been no contamination to the parts touching the patient. General precautions must be observed in order to avoid cross-contamination of disease from one patient to another. If ear cushions or ear tips are contaminated, it is strongly recommended to remove them from the transducer before they are cleaned. Water should be used for frequent cleaning, but in the case of severe contamination it may be necessary to use a disinfectant. The use of organic solvents and aromatic oils must be avoided.
9. Great care should be exercised when handling transducers, as rough handling, for example dropping onto a hard surface may break or damage the parts.


1.4 Concerning Repair

Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:


- assembly operations, extensions, readjustments, modifications or repairs are carried out by authorized persons,
- a 1 year service interval is maintained
- the electrical installation of the relevant room complies with the appropriate requirements, and
- the equipment is used by authorized personnel in accordance with the documentation supplied by Interacoustics.

It is important that the customer (agent) fills out the RETURN REPORT every time a problem arises and sends it to Interacoustics, Drejervaenget 8, DK-5610 Assens, Denmark. This should also be done every time an instrument is returned to Interacoustics. (This of course also applies in the unlikely worst case scenario of death or serious injury to a patient or user).

1.5 Titan for Hand Held Use (IMP440)


1.6 Titan for Hand Held Use (DPOAE440)


Manual Start (DPOAE440)

1. Insert probe
2. Select ear
3. Pres **Start**
4. View recording & wait for **Done!**
5. Insert probe in other ear
6. Select other ear
7. View recording & wait for **Done!**

Select client & Save

1.7 Titan for Hand Held Use (ABRIS440)

Before Start

1. Prepare patient
2. Place electrodes
 - White = vertex (forehead)
 - Red = ground (cheek/shoulder)
 - Blue = nape

3. Place transducer

4. Check electrode impedances

14.1 KΩ 4.2 KΩ 43.8 KΩ

Green = Acceptable
Red = Poor

5. Select ear (for probe)

The binaural ear indicator will appear if a binaural transducer and nape montage is used.


6. Press **Start**

7. View recording & wait for **Done!**

8. If using probe or mastoid montage with binaural transducer select other ear and/or place transducer.

9. Press **Start**

10. View recording & wait for **Done!**

Edit before saving

1.8

1.8 Titan for PC-Controlled Use (IMP440 & DPOAE440)

1. Place Titan in its cradle. Use extension cord with shoulder box

2. Switch Titan on by selecting ear


4. Start OtoAccess™ or Noah3, select patient, open Titan Suite

3. Attach the contra head phone or insert to the shoulder box if contra reflex testing is to be performed (IMP only).

15. Print data

16. Save data

5. Select the Module tab for the relevant test

6. Select protocol

7. Select ear

8. Place transducer/s in the ear/s

9. Press Start

10./14. View recording and wait until done


11. Select other ear

12. Place transducer/s in the ear/s

13. Press Start


Swap combined / single display (IMP)




Toggle monaural/ binaural view (DPOAE)

Review data

Check if need to retest

1.9 Titan for Pc-Controlled Use (ABRIS440)

Note: When using a binaural transducer (inserts or headphones) and the nape montage, both ears can be tested simultaneously.


1.10 Creating a new Protocol

Creating a new protocol

1. Go to the required **Module** tab.
2. Select **Menu | Setup | Protocol Setup**.
3. Press **New**.
4. Type new name.
5. Select test to include in the protocol.
6. Press **Add** to include selected tests.
7. Reorder tests if necessary.
8. Press **Settings** and thoroughly go through the different settings of all tests. Do not forget advanced settings.
9. Optionally, personalize displays, press **Display wizard**.
10. Optionally, select print template, press **Print wizard**.
11. Save protocol, press **Ok**.


Optionally, transfer protocol from Titan Suite to Titan (see Quick guide: Titan data flow).

Note: Protocol setup screens may differ slightly between modules


Using Display Wizard

1. Press right mouse button on an empty space to add more displays or a combined display.
2. Select how to use the combined display.
3. Select the display you want to modify.
4. Select during which tests the current display is used.
5. Press right mouse button on empty space to add a certain graph.
6. Use right mouse click on graph to change settings.
7. Drag the graph over the screen and onto the drop indicators to give it its appropriate place in the display.
8. Press **OK** to save when done.



Using Print Wizard

1. Select an existing template or
2. Press **New template / Edit template**.
3. Add pages to the print template.
4. Select which page you want to modify.
5. Drag items from the lists onto the page
6. Resize and place items on the page
7. Use right mouse click on items to change settings.
8. Press **Save** or **Save as...** when done.
9. Select the created template.
10. Press **Select** in order to automatically use the template for this protocol.


1.11 Titan Data Flow

OtoAccess™ or Noah3 database

Transfer protocols


1. Select Menu | Setup | Hardware protocols at Main tab.
2. Select protocols from Protocol Library
3. Press Add
4. Arrange order on hardware
5. Press OK

Download sessions to database

1. Press **Session download** at Main tab
 2. Select which sessions need **Copy to PC**
 3. Press **Download**
- Known (uploaded) clients will not prompt
4. Sessions from unknown clients require a search in your database...
 5. ...or press **Enter new client**
 6. Press **Save**

Note: Test screen displays will differ between modules

Upload client names and details

1. Press **Client upload** at Main tab
2. Search patients in your database
3. Press **Add**
4. When all required patients are found and added, press **Download to hardware**

1.12 Titan - Operation Panel Dictionary

Position	Symbol	Function
1 - 3	F1 – F3	The F1 – F3 keys allow selection of various options indicated on the Titan display screen above each individual function key, e.g., “Print”, “Save” and “Delete”
4 + 5	F4 – F5	The R and L buttons allow ear selection. Hold down F4 (R) and F5 (L) at the same time to turn off the Titan.
6 + 7	F6 – F7	The Up and Down buttons (F6 and F7) navigate through lists, menus etc. These buttons also allow you to scroll through the different tests which are indicated in the selected screen.

1.13 Technical Specifications

Medical CE-mark:	The CE-mark indicates that Interacoustics A/S meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. Approval of the quality system is made by TÜV – identification no0123	
Standards:	Safety:	IEC 60601-1, Internally powered, Type B and BF applied parts
	EMC:	IEC 60601-1-2
	Impedance:	IEC 60645-5/ANSI S3.39, Type 1
	Test Signal:	IEC 60645-1/ANSI S3.6 , IEC 60645-3
	OAE:	IEC 60645-6 2009, Type 2 Otoacoustic emissions
	ABR:	IEC 60645-7 2009, Type 2
Cradle	Safety:	IEC 60601-1, Class II
	Power	Astrodyne ASA30M-0301
	Mains voltages and frequencies:	100 – 240 VAC, 47 – 63 Hz
	Consumption:	0.8 – 0.4 A
Operation environment:	Temperature:	15 – 35 °C
	Rel. Humidity:	30 – 90%
Storing/handling:	Temperatures below 0°C and above 50°C may cause permanent damage to the instrument and its accessories.	
Transport:	Temperature:	-20 – 50 °C
	Rel. Humidity:	10 – 95%

Impedance Measuring System

Probe tone:	Frequency:	226 Hz, 678 Hz, 800 Hz, 1000 Hz
	Level:	69 dB HL with AGC, assuring constant level at different ear canal volumes.
Air pressure:	Control:	Automatic.
	Indicator:	Measured value is displayed on the graphical display.
	Range:	-600 to +300 daPa.
	Pressure limitation:	-750 daPa and +550 daPa.
	Pressure change rate:	Minimum, medium, maximum or automatic with minimum speed at compliance peak. Selectable in the setup.
Compliance:	Range:	0.1 to 8.0 ml at 226 Hz probe tone (Ear volume: 0.1 to 8.0 ml) and 0.1 to 15 mmho at 678, 800 and 1000 Hz probe tone.
Test types:	Tympanometry	Automatic, where the start and stop pressure can be user-programmed in the setup function. Manual control of all functions.
	Eustachian tube function 1 – Non perforated eardrum	Williams test
	Eustachian tube function 2 – Perforated eardrum	Toynbee test
	Eustachian tube function 3 – Patulous Eustachian tube	Continuous sensitive impedance measurement for 30 up to 150 s.
	Indicators:	Graphical display Compliance is indicated as ml and pressure as daPa. In PC controlled mode admittance, susceptance and conductance can be printed. Stimulus level is indicated as dB Hearing Level.
Memory:	Tympanometry:	1 curve per ear per tympanometry test. 3 curves per ear per Eustachian tube function test. And theoretically an infinite number of tests per protocol.

There is no deviation between static and dynamic mode.

Reflex Functions

Signal sources:	Tone - Contra, Reflex:	250, 500, 1000, 2000, 3000, 4000, 6000, 8000 Hz.
	Tone - Ipsi, Reflex:	500, 1000, 2000, 3000, 4000 Hz.
	NB noise - Contra, Reflex:	250, 500, 1000, 2000, 3000, 4000, 6000, 8000 Hz.
	NB noise - Ipsi, Reflex:	1000, 2000, 3000, 4000 Hz.
	Noise - Contra, Reflex:	Wide Band, High Pass, Low Pass.
Outputs:	Noise - Ipsi, Reflex:	Wide Band, High Pass, Low Pass.
	Contra Earphone:	TDH39 earphone, DD45 earphone, CIR55 insert and/or EARTone 3A insert for Reflex measurements.
	Ipsi Earphone:	Probe earphone incorporated in the probe system for Reflex measurements.
	Air:	Connection of the air system to the probe.
	Test types:	
	Manual Reflex	Manual control of all functions.
	Automated Reflex	Automatic reflexes: - Single intensities - Reflex growth
	Reflex Decay	Automatic, 10 dB above threshold and manually controlled with stimulus durations of 10 to 30 s.
	Reflex latency	Automated, first 300 ms from stimulus start.

ABR Infant Screening

Preamplifier	One Channel	3 electrodes. 50 cm Switchable: Software will automatically switch mastoid and ground if mastoid montage is used. So the user does not need to replace electrode during testing.
Electrical Impedance measurement	Gain:	64 dB
	Frequency response:	0,5 - 5000 Hz
	Noise:	<25 nV/ $\sqrt{\text{Hz}}$
	CMR Ratio:	>90 dB.
	Max input offset voltage:	2.5 V
Stimulus	Input impedance:	10 M Ω / 170 pF
	Power from main unit:	Isolated power supply
	Measurement frequency:	33 Hz
	Waveform:	Rectangular
	Measurement current:	11.25 μA
Recording	Range:	0.5 k Ω – 25 k Ω $\pm$ 10 %
	Stimuli	Click range (177-11313 Hz) CE-Chirp® range (177 – 11313 Hz) HiLo CE-Chirp® range (Lo 177- 1414 Hz) & (Hi 1414 – 11313 Hz)
	Stimulus rate:	90 Hz
	Transducers:	EarTone 3A insert phone auto calibrated & detection TDH 39 head phone auto calibrated & detection IOW Probe auto calibrated & detection
	Channels:	2
Display	Level:	30 dB nHL, 35 dB nHL, 40 dB nHL
	Bandwidth:	22.05 kHz
	Analysis time:	1-10 min or Residual noise 5-80 nV
	A/D resolution:	24 bit
	Artifact reject system:	Rejection level (Peak, Min RMS, Max RMS) Clipping (Saturation)
Algorithmic Sensitivity	CE-Chirp®	99.9%
Algorithmic Specificity	CE-Chirp®	> 96%

DPOAE

Stimulus	Frequency range:	500 to 10000 Hz
	Frequency step:	25 Hz
	Level:	30 to 80 dB SPL (75 dB SPL for 6kHz and 65 dB SPL for 8kHz to 10kHz)
	Level Step:	1 dB
Recording	Transducer:	IOW Probe auto detection, auto calibrated
	Analysis time:	Minimum 2 seconds to unlimited time
	A/D Resolution:	24 bit, 5.38 Hz resolution
	Artifact rejection system:	-30 to +30 dB SPL or off
	Stimulus tolerance:	Adjustable between 1 and 10 dB
	SNR criteria:	Adjustable between 3 and 25 dB
	Probe check window:	256 points frequency response of the ear canal due to a click stimulus.
Display	DP-response window:	4096 points frequency response
	General display gain:	Applicable during testing
	Display	Stimulus level and type, Bar and Graph view
Probe specifications	Titan IOW probe:	IMP, DPOAE and ABRIS capable Replaceable probe tip
Other		
Test Pressure:		Ambient pressure.

General

PC control:	USB:	Input/output for computer communication. Titan can be fully operated from a PC. The measurements can then be followed on the PC screen. Data can be sent to and saved on the PC and stored in OtoAccess™. See separate section in Service Manual for programming details.
Memory:		Theoretically, an infinite amount of test results can be stored on the PC. The Titan hand held unit is delivered with a 1 Gb memory card, enough for storing more than a quarter of a million tests.
Thermal printer (Optional):	Type:	Thermal (Bluetooth) printer with recording paper in rolls. Print on command through Bluetooth communication and through serial RS-232.
	Paper width:	57.5 ± 0.5 mm on thermal printer
	Printing time:	Printing time depends on the size of the used protocol. For 2 tympanograms and 8 reflexes the thermal printer uses approximately 6s.
Dimensions		6x6x28cm/2.4x2.4x11inch
Titan Weight		360g/0.8lbs
ABR/OAE/IMP shoulderbox weight		120 g
ABR/OAE/IMP shoulderbox Dimensions		102, 68, 26 mm.
OAE/IMP shoulderbox weight		64 g
OAE/IMP shoulderbox Dimensions		65, 35, 18 mm.

Table 1: Frequencies and Intensity Ranges for IMP440

Titan Maximums IMP										
Center Freq.	TDH39		CIR55		EARTone 3A		IOW IPSI		DD45	
	Reading		Reading		Reading		Reading		Reading	
	Tone	NB	Tone	NB	Tone	NB	Tone	NB	Tone	NB
	[Hz]	[dB HL]	[dB HL]	[dB HL]	[dB HL]	[dB HL]	[dB HL]	[dB HL]	[dB HL]	[dB HL]
125	80	65	85	70	100	85	70	60	80	65
250	100	85	100	85	110	100	85	75	100	85
500	120	100	110	100	115	105	100	85	115	100
750	120	105	110	105	120	110	100	85	120	105
1000	120	105	115	105	120	110	105	90	120	105
1500	120	105	115	105	120	110	110	90	115	100
2000	120	105	115	105	120	110	105	90	115	100
3000	120	105	115	105	120	110	95	90	125	105
4000	120	105	110	100	120	105	100	85	115	105
6000	120	100	95	95	105	100	85	80	110	90
8000	105	95	80	85	90	85	80	75	105	95
10000										
WB	-	120	-	120	-	120	-	105	-	125
LP	-	120	-	120	-	120	-	110	-	120
HP	-	120	-	120	-	120	-	105	-	130

Table 2: Frequencies and Intensity Ranges for DPOAE440

Titan Maximums DPOAE		
Center Freq.	IOW IPSI	IOW ch2
	Reading	Reading
	Tone	Tone
[Hz]	[dB SPL]	[dB SPL]
500	80	80
750	80	80
1000	80	80
1500	80	80
2000	80	80
3000	80	80
4000	80	80
6000	75	75
8000	65	65
10000	65	65

ABRIS maximum levels are not listed as these are limited to 30, 35 & 40 dBnHL for all transducers.

Specification of Input/Output Connections

Inputs	Connector type	Electrical properties	
Patient response	Jack, 3,5mm 4-pole	Handheld switch: Pin 1: GND Pin 2: Signal Pin 3: Future use I/O Pin 4: Future use I/O	3V through 10K Ω is forced to ground when activated
Outputs:			
Phones, Left/ Right	Jack, 3,5mm 4-pole	Voltage: Min. load impedance: Pin 1: CH1 GND Pin 2: CH1 OUT (left) Pin 3: CH2 OUT (right) Pin 4: CH1 GND	Up to 3V rms. by 10 Ω load 8 Ω Pin 3:
Phones, Contralateral	Jack, 3,5mm 4-pole	Voltage: Min. load impedance: Pin 1: CH1 GND Pin 2: CH1 OUT (left) Pin 3: CH2 OUT (right) Pin 4: CH1 GND	Up to 3V rms. by 10 Ω load 8 Ω
Transducer	IA proprietary, 12-pole	Pin 1: Pin 2: Pin 3: Pin 4: Pin 5: Pin 6: Pin 7: Pin 8: Pin 9: Pin 10: Pin 11: Pin 12:	CH1 out CH1 GND DGND GND A / GND Microphone Microphone – input / Analog balanced in Microphone + input / Analog balanced in Power supply +3/+5V CH2 out CH2 GND I2C CLK I2C DATA I2C Interrupt
Data I/O: USB	USB type"B"	USB port for communication	See appendix A in service manual for detailed information

Calibration Properties

Calibrated Transducers:	Contralateral Earphone:	Telephonics TDH39 with a static force of 4.5N $\pm$ 0.5N and/or EARTone 3A and/or CIR55 insert phone
	Probe system:	Ipsilateral Earphone: is integrated in the probe system Probe frequency transmitter and receiver and pressure transducer is integrated in the probe system
Accuracy:	General	Generally the instrument is made and calibrated to be within and better than the tolerances required in the specified standards:
	Reflex Frequencies:	$\pm$ 1%
	Contralateral Reflex and Audiometer Tone Levels:	$\pm$ 3 dB for 250 to 4000Hz and $\pm$ 5 dB for 6000 to 8000Hz
	Ipsilateral Reflex Tone Levels:	$\pm$ 5 dB for 500 to 2000Hz and +5/-10 dB for 3000 to 4000Hz
	DPOAE Levels:	$\pm$ 1.5 dB for 1000 to 4000Hz and $\pm$ 3 dB outside range
	Pressure measurement :	$\pm$ 5% or $\pm$ 10 daPa, whichever is greater
	Compliance measurement:	$\pm$ 5% or $\pm$ 0.1 ml, whichever is greater

Impedance Calibration Properties

Probe tone	Frequencies:	226 Hz $\pm$ 1%, 678 Hz $\pm$ 1%, 800 Hz $\pm$ 1%, 1000 Hz $\pm$ 1%
	Level:	85 dB SPL $\pm$ 1.5 dB measured in an IEC 60318-5 acoustic coupler. The level is constant for all volumes in the measurement range.
	Distortion:	Max 1% THD
	Range:	0.1 to 8.0 ml

Compliance	Temperature dependence:	-0.003 ml/°C
	Pressure dependence:	-0.00020 ml/daPa
	Reflex sensitivity:	0.001 ml is the lowest detectable volume change
	Temporal reflex characteristics:	Initial latency = 35 ms (±5 ms)
		Rise time = 45 ms (±5 ms)
		Terminal latency = 35 ms (±5 ms)
		Fall time = 45 ms (±5 ms)
		Overshoot = max. 1%
Undershoot = max 1%		
Pressure	Range:	Values between -600 to +300 daPa can be selected in the setup.
	Safety limits:	-750 daPa and +550 daPa, ±50 daPa

Reflex Calibration Standards and Spectral Properties:

General Specifications for stimulus and audiometer signals are made to follow IEC 60645-5

Contralateral Earphone	Pure tone:	ISO 389-1 for TDH39 and ISO 389-2 for CIR 22.
	Wide Band noise (WB):	Interacoustics Standard
	– Spectral properties:	As “Broad band noise” specified in IEC 60645-5, but with 500 Hz as lower cut-off frequency.
	Low Pass noise (LP):	Interacoustics Standard
	– Spectral properties:	Uniform from 500 Hz to 1600 Hz, ± 5 dB re. 1000 Hz level
Ipsilateral Earphone	High Pass noise (HP):	Interacoustics Standard
	– Spectral properties:	Uniform from 1600 Hz to 10KHz, ± 5 dB re. 1000 Hz level
	Pure tone:	Interacoustics Standard.
	Wide Band noise (WB):	Interacoustics Standard
	– Spectral properties:	As “Broad band noise” specified in IEC 60645-5, but with 500 Hz as lower cut-off frequency.
	Low Pass noise (LP):	Interacoustics Standard
	– Spectral properties:	Uniform from 500 Hz to 1600 Hz, ± 10 dB re. 1000 Hz level
	High Pass noise (HP):	Interacoustics Standard
	– Spectral properties:	Uniform from 1600 Hz to 4000 Hz, ± 10 dB re. 1000 Hz level
	General about levels:	The actual sound pressure level at the eardrum will depend on the volume of the ear. See Table 2 for details.

The risk of artifacts at higher stimulus levels in reflex measurements are minor and will not activate the reflex detection system

Titan Instruction for Use - English

Date: 2011-06-06

Page 20/21

Table 3: Reference Values for Stimulus Calibration

Freq.	Reference values for stimulus calibration [dB re. 20 µPa]								Variation of Ipsi stimulus levels for different volumes of the ear canal Relative to the calibration performed on an IEC 126 coupler [dB]		Sound attenuation values for TDH39 earphones using MX41/AR or PN51 cushion [dB]
		ISO 389-1 (Interacoustics Standard)	ISO 389-2 (Interacoustics Standard)	ISO 382-2 (Interacoustics Standard)	Interacoustics Standard	Interacoustics Standard	Interacoustics Standard	ISO 389-4 (ISO 8798)	0.5 ml	1 ml	
[Hz]		TDH39	EARTone 3A / ABR	CIR55	DD45	IOW Probe	IOW Probe NB	NB Stimulus Correction Values (except IOW probe)			
RETSP1	125	45	26	26	47.5	41	43.5	4			3
	250	25.5	14	14	27	24.5	26.5	4			5
	500	11.5	5.5	5.5	13	9.5	17	4	9.7	5.3	7
	1000	7	0	0	6	6.5	10.5	6	9.7	5.3	15
	1500	6.5	2	2	8	5	12	6			21 (1600 Hz)
	2000	9	3	3	8	12	11	6	11.7	3.9	26
	3000	10	3.5	3.5	8	11	11	6	-0.8	-0.5	31 (3150 Hz)
	4000	9.5	5.5	5.5	9	3.5	8	5	-1.6	-0.8	32
	6000	15.5	2	2	20.5	3	5.5	5			26 (6300 Hz)
	8000	13	0	0	12	-5	-0.5	5			24
	WB	-8	-5	-5	-8	-5			7.5	3.2	
	LP	-6	-7	-7	-6	-7			8.0	3.6	
	HP	-10	-8	-8	-10	-8			3.9	1.4	
neRETSP1	CE-Chirp®	27.5	31.5		26	32					
	CE-Low Chirp®	26.5	26.5		25.5	27.5					
	CE-High Chirp®	28	31		28	32					
	Click	30.5	35		32.5	28					

**All figures in bold are Interacoustics Standard values.*

Coupler Types used by Calibration

IMP:

TDH39 is calibrated using a 6cc acoustic coupler made in accordance to IEC 60318-3, Ipsilateral earphone and probe tone are calibrated using a 2cc acoustic coupler made in accordance to IEC 60318-5

DPOAE:

Probe stimuli L1 and L2 are calibrated in SPL values using an ear simulator coupler made in accordance to IEC 60318-4.

ABRIS:

Probe and insert stimuli are calibrated in SPL values using an ear simulator coupler made in accordance to IEC 60318-4. Headphones (TDH39 and DD45) stimuli are calibrated in SPL values using an artificial ear coupler according to IEC 60318-1.

General about Specifications

Interacoustics continuously strives to improve its products and their performance. Therefore the specifications can be subject to change without notice.

The performance and specifications of the instrument can only be guaranteed if it is subject to a technical maintenance at least once a year. This should be carried out by a workshop authorized by Interacoustics.

Interacoustics puts diagrams and service manuals at the disposal of authorized service companies.

Enquiries about representatives and products may be sent to:

Interacoustics A/S	Phone:	+45 63713555
Drejervaenget 8	Fax:	+45 63713522
DK 5610 Assens	e-mail:	info@interacoustics.com
Denmark	http:	www.interacoustics.com

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The diagram illustrates the workflow for using the Titan device. It includes screenshots of the 'Select Protocol', 'Test: Out of ear', and 'Done! DP-Gram' screens. Arrows indicate the sequence of actions: selecting a protocol, starting the test, viewing results, and printing. A list of 7 steps for manual start is provided, along with a 'Select Client & Save' screen showing a list of clients and an 'Edit New' screen for editing client information.

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